



GYN-3000E - Electric Table for Gynaecology and Obstetrics

GYN-3000E



GYN-3000E Electric Table for Gynaecology and Obstetrics

Tabletop Length/Width	1780/600mm
Tabletop Up/Down	880/680mm
Tabletop Trendelenburg/Reverse Trendelenburg	8°/20°
Back Plate Up/Down	45°/7°
Voltage/Frequency/Capacity	220V/50Hz/1.0kw
Emergency Battery	Optional

Introduction

GYN-3000E Is the most latest electric table that can be used for Obstetrics birthing, Gynecological operation and examination

The tabletop elevation, trendelenburg/reverse trendelenburg, back plate all controlled by both micro touch hand remote control and foot switch

Imported electric motors enable flexible, smooth movements, low noise, high reliability

All base and column cover are made of 304 medical grade stainless steel

Seamless mattress can be colorfully customized, easy to clean and sterilize