

Infant Phototherapy PH-10

Function features

- Equal light distribution
- High intensity
- Quiet, no noise
- Five blue fluorescent lights
- The angle of head and height are adjustable
- Steel support pole, four castors with brake
- Integral time record therapy time



Technical data

Specification	PH-10	
Power Supply	AC 110/220V 60/50Hz	
Power input	120VA	
Wavelength (blue light)	400nm ~ 550nm	
Total Irradiance for bilirubin	1300 μ W/cm ²	
Height Adjust Range	1100 ~ 1600mm	
Pitching angle of phototherapy head	$\pm 90^\circ$	
Time Accuracy	1min/12h	
Integral time range	0h ~ 9999h 59min	
Operating Condition	Temperature	18 $^\circ$ C ~ 30 $^\circ$ C
	Relative Humidity	10% ~ 85%RH
	Atmospheric pressure	700hpa ~ 1060hpa
Transport and storage	Ambient temperature	-10 $^\circ$ C ~ +55 $^\circ$ C
	Relative humidity	$\leq 95\%$
	Atmospheric pressure	500hpa ~ 1060hpa