

Vital Signs Patient Monitor SCOPE700

Specifications



Oxford Medical Group Inc.

869 E. Schaumburg Rd. Suite 194
Schaumburg, IL 60194 - USA
Phone; +630-847+0744
E-mail: sales@oxford.com

Features:

- ◆ 7" high resolution color TFT LCD display.
- ◆ Real time S-T segment analysis, pace-maker detection.
- ◆ Indication of Incorrect manipulation and malfunction analysis.
- ◆ Patient information input management function
- ◆ Multi-lead ECG waveforms display in phase.
- ◆ Large volume of tabular and graphic trends information storage and easy to recall.
- ◆ Capture dynamic waveforms.
- ◆ Efficient resistance to interference of defibrillator and HF knife.
- ◆ 4 hours working capacity with built-in rechargeable lithium battery.
- ◆ Unique accessory management cabinet.
- ◆ Networking capacity.

Technical Parameters:

- ◆ 4 types configuration is available according to users request:
A type: SpO2, NIBP, HR
B type: SpO2, NIBP, HR, ECG, RESP
C type: SpO2, NIBP, HR, TEMP
D type: SpO2, NIBP, HR, ECG, RESP, TEMP
- ◆ Option: Printer, ETCO2, IBP, wall mount, trolley
- ◆ Machine Size: 17.5*14*30cm N.W: 3.3 kg

Packing Info:

- main unit
- Disposable ECG Electrodes 5 pcs
- Reusable SpO2 sensor 1 PC
- ECG Cable and Leadwire 1 PC
- NIBP cuff and tube 1 PC
- Built-in lithium battery 1 PC
- Rechargeable lithium battery 1 PC
- Power cable 1 PC
- Battery 1 PC