

# SMART-1200A

12-channel Electrocardiograph



- Twelve Channel ECG Machine
- 8.4 inch (800x600) high-resolution touch screen
- Sleep/wake up mode for saving energy and extending LCD and battery life
- Last 120 seconds waveforms reviewable and recordable
- Internal storage for up to 300 digital ECG records
- External printer support and USB socket for data transmission
- 122 kinds of arrhythmia analysis
- Pace-maker detection
- Alphanumeric keyboard, single button operation for sampling, analyzing, printing, storage and transfer



**Oxford Medical Group Inc.**

869 E. Schaumburg Rd. Suite 194

Schaumburg, IL 60194 - USA

Phone: +630-847-0744

E-mail: [sales@ocfford.com](mailto:sales@ocfford.com)

