

SMART-1200B


12-channel Electrocardiograph



Features

- 5.7 inch TFT color screen
- Lead status indication
- Simple & clear keyboard layout, back light support
- Internal memory of 300 ECG records
- High accurate digital filter (EMG filter, AC filter, Drift filter, lowpass filter)
- 122 kinds of diagnosis report
- 120 seconds ECG waveforms review and print
- Support 210mm rolling, Z-folded paper
- USB port for external printer connection & data transmission
- PC-ECG management system
- Working mode: Auto, Manual, Rhythm

Interface



Various printing mode



Printing mode: 12X1, 6X2+1R, 3X4+1R, auto analysis



SMART-1200B

12-channel Electrocardiograph



Oxford Medical Group Inc.

869 E. Schaumburg Rd. Suite 194
Schaumburg, IL 60194
USA

Phone: +630-847-0744 / E-mail: sales@oxford.com



SMART-1200B

12-channel Electrocardiograph



Hidden handle




Alphanumeric keyboard;
Back light for button



Memory of 300 ECG records;
USB for data transmission;
PC-ECG management system

Clear professional report



5.7 inch TFT color screen,
easy for observation



Internal aluminum alloy rack ,
shockproof, Anti-interference
design



High quality printing
system , print is clear and
smooth



Internal modular design,
easy for maintenance

