

SONO-OX50 Touch

- ❖ 12.1" folding-up touch screen, 0~90°C multi-angle observation.
- ❖ Parameter: ECG/HR, RESP, NIBP, SpO2, PR, TEMP, TOCO, FHR, FM.
- ❖ Various interfaces: main, fetal monitoring, fetal monitoring cascade, big font, multi-parameter, trend graphic display.
- ❖ Manual/automatic fetal movement measurement, automatic fetal status scoring system.
- ❖ 112mm thermal printer, real-time/timer printing for fetal waveforms and parameters.



SONO-OX50 PLUS

- ❖ 6.5" color TFT touch screen, convenient for observation & movement. Built-in 152mm width printer, complied with international standard.
- ❖ Automatic detection transducers, flexible connection with sensor, self-optimized interfaces.
- ❖ 12-crystal wide pulse wave waterproof transducer, suitable for underwater delivery.
- ❖ Advanced patient files management system.
- ❖ TEMPOMTM timing monitoring function, avoid overtime monitoring. Intelligent alarm management, automatic-identification of alarm level. Self-adjust proper alarm time to reduce false alarms.



SONO-OX50

- ❖ 5.7" LCD screen.
- ❖ Built-in high-capacity Li-ion battery, above 4 hours working time.
- ❖ 112mm wide printer. Special extractable paper slot, easy for paper replacement and observation.
- ❖ Support 30mm/Min, 20mm/Min, 10mm/Min real-time printing, and 25mm/Sec high-speed fetal waveform review printing.
- ❖ Support timing printing.
- ❖ 3-level audio/visual alarms.

Standard Configuration

- ❖ FHR, TOCO, FM

Optional Configuration

- ❖ Twin, Fetal stimulator

